








# OGUK

MEDICAL - ORTHOPAEDIC - COMPRESSION - DIABETIC



**OGUK**

**MEDICAL STOCKINGS**

## PRODUCT QUALITY

# 100% Made in Italy

The standard features of the textile production process, which includes the employment of selected yarns and effective production methods resulting from long experience in this field, guarantee:

- Elasticity
- 4-way stretch
- Fabric elongation
- Durability
- Comfort and wearability
- Shape-memory effect

A highly-skilled and well-trained staff carefully checks the product during the whole production process:

- Just after the production phase (weaving)
- During and after the product assembly phase (sewing)
- During the vapor ironing
- After the dyeing process
- During the packaging process
- Before dispatching the product

## CERTIFIED QUALITY

# Product and System certificates

- Certified graduated compression in compliance with the officially defined and recognized medical standards, both Italian and international: CE Medical Device, FDA, Oeko-Tex, ISO 9001 and ISO 13485.
- Guaranteed administration of mmHg at ankle.
- Proven efficacy, effective, repetitive, permanent and extended results.
- Anatomical shapes obtained during the production process, which prevent the product from being too much tight on knees and thighs.



## VENOUS INSUFFICIENCY

# Therapeutic dosage program of mmHg

The appropriate and personalized administration of mmHg is both the basic treatment for venous and lymphatic system diseases and an adjuvant to other types of treatment and prophylaxis.

These dosages are measured in mmHg (i.e. value/power of the compression at ankle) and turn into an external compression of veins.

As a result:

- Enlarged veins are compressed until their normal tone and diameter are restored.
- Venous valves are restored to their normal functions and prevent the blood from flowing backwards.
- Venous blood flows more rapidly back to the heart.
- Leg swelling, heaviness and pain decrease.
- Venous blood pressure drops.
- The risk of blood clot formation decreases.

The dosage in mmHg varies according to the clinical manifestation of venous insufficiency and shall be increased in relation to the progressive appearance of symptoms - from swollen ankles to skin blemishes, telangiectasias, superficial or deep thrombophlebitis, and ulcers.

The progression of venous insufficiency is mainly due to lifestyle, sedentariness, dietetic errors, age, and heredity.

**OGUK** MEDICALE products come in a variety of dosages (i.e. from 15 to 46 mmHg), lengths, and shapes. The application of the most suitable product makes a crucial contribution toward the treatment for venous insufficiency and the slowdown in its progression without any contraindications. Furthermore, it helps to treat venous and lymphatic system diseases.

IS ESSENTIAL TO CHOOSE THE RIGHT SIZE of the garments that patients should wear in order to guarantee the correct administration of mmHg. The size shall be chosen according to the ankle circumference and not to the clothing size.

IF THE SIZE IS WRONG, THE DOSAGE WILL BE UNSUITABLE; as a consequence, the result will not be satisfactory.

## THE **OGUK** MEDICALE SOLUTION

# Compression Categories

The therapeutic graduated compression stockings are subdivided into compression Classes according to the experimental UNI ENV 12718 European standard depending on the compression strength at the ankles. The therapeutic graduated compression stockings are normally prescribed by specialized physicians who

recommend the most suitable garment and compression for the pathology to be treated.

The 3 different compression classes of the therapeutic **OGUK** Medicale stockings are generally prescribed for the following venous and lower extremity disease:

Clinical Picture	Recommended Compression	Solution
Prophylaxis of venous diseases in risk patients, tired and heavy sensation in the legs, appearance of varicose veins during pregnancy, thrombosis and embolism.	<b>Medium compression</b> CLASS 1 - K1 - 15-21 mmHg* Anti-Embolism AE20 - 18-23 mmHg*	CLASSIC SOFT COTTON ANTI-EMBOLISM
Slight chronic venous insufficiency, thrombophlebitis (superficial thrombosis) with and without oedema, post-sclerotherapy, after selective out-patient treatment of varices and stripping, varicose veins during pregnancy, prophylaxis of thrombosis and embolism.	<b>Moderate compression</b> CLASS 2 - K2 - 23-32 mmHg* Anti-Embolism AE30 - 25-32 mmHg*	CLASSIC SOFT COTTON ANTI-EMBOLISM
Advanced phases of chronic venous insufficiency, after recovery from crural ulcer (prophylaxis of relapse), post-traumatic oedema, and reversible lymphedemas.	<b>Strong compression</b> CLASS 3 - K3 - 34-46 mmHg*	CLASSIC

\*Mentioned compression refers to the ankle and decreases gradually to the groin, 70% on the calf and 40% on the thigh, guaranteeing a right and good blood circulation.

# CLASSIC

Over 20% of the adult population suffer from venous insufficiency and show increasingly severe symptoms, which may include:

- TIRED LEGS AND SWOLLEN ANKLES
- SKIN BLEMISHES
- OEDEMA
- SUPERFICIAL AND DEEP THROMBOPHLEBITIS
- SKIN LESIONS.

In all these cases, the **OGUK** MEDICALE CLASSIC line is not only a basic treatment for venous hypertension but also an adjuvant to other types of treatment and prophylaxis.



**OGUK** MEDICALE CLASSIC products come in a variety of lengths and types – normal or short length, closed or open toe – in order to provide patients with

comfort and to guarantee an excellent result. In addition to standard sizes, the manufacturing of made-to-measure articles is also possible.

It is possible to choose the most suitable dosage of mmHg for every patient. Three types of products are available and they all have a high percentage of elastic yarns and exert graduated compression:

Compression	AD	AG	AGH	AG TL/R	AT	ATM	Sizes*	Length		Toe		Colors	Composition
CLASS 1 - K1 15-21 mmHg	✓	✓	✓	✓	✓	✓	From 1 to 5	Normal	Short	Closed	Open	Beige	75% PA 25% EA
CLASS 2 - K2 23-32 mmHg	✓	✓	✓	✓	✓	✓	From 1 to 5	Normal	Short	Closed	Open	Beige	70% PA 30% EA
CLASS 3 - K3 34-46 mmHg	✓	✓	✓		✓		From 1 to 5	Normal		Closed	Open	Beige	70% PA 30% EA

\*Size Table on the last page



# ANTI-EMBOLISM

The blood flows back to the heart thanks to four “pumps”: the foot pump, the calf muscle pump, the thoracic pump, and the cardiac pump. When one or more of these pumps do not work properly, venous blood becomes stagnant in the legs and may cause a considerable reduction in blood circulation speed, tendency to venous thrombosis, and possible thromboembolism. Furthermore, if venous circulation is already damaged due to varicose veins, the risk of thrombosis considerably increases. Indeed, thromboembolism is the third cause of cardiovascular death, after heart attack and stroke.



**OGUK** MEDICALE ANTI-EMBOLISM prevents the risk of thromboembolism, if you administrate 18-23 mmHg or 25-32 mmHg beforehand. Then, it is possible to continue with a personalized application of **OGUK** MEDICALE CLASSIC, SOFT, or COTTON products.

**OGUK** MEDICALE ANTI-EMBOLISM products have an excellent effectiveness because they exert graduated compression on legs and are not too much tight over the ankle.

Compression	AD	AGH	AG TL/R	Sizes*	Length	Toe	Colors			Compo- sition	Packaging	
AE20 STANDARD 18-23 mmHg	✓	✓		From S to XXL	Normal	Open	White	Beige	Black	75% PA 25% EA	Polybag	Box
AE20 PREMIUM 18-23 mmHg	✓	✓	✓	From S to XXL	Normal	Inspect- able	White	Beige	Black	75% PA 25% EA	Polybag	Box
AE30 STANDARD 25-32 mmHg	✓	✓		From S to XXL	Normal	Open	White	Beige	Black	70% PA 30% EA	Polybag	Box



Dispositivo Medico  
Medical Device



14CX00025 CENTROCOT  
Tested for harmful substances.  
www.oeko-tex.com/standard100

# CLASSIC ARM SLEEVE

The line of compression arm sleeves provides patients with all the compression and qualitative properties of the **OGUK** MEDICALE CLASSIC line for lower limbs.

The **OGUK** compression arm sleeves are generally prescribed for the following diseases:


- LYMPHOEDEMA AT AN EARLY STAGE
- SWELLING AFTER BREAST CANCER SURGERY
- TREATMENT FOR BURNS AND WIDE SCARS
- SWELLING AFTER INJURIES OR SURGERY.




Products are available in the following versions:

- Basic compression arm sleeves with or without hands
- Compression arm sleeves with bra strap hooks, with or without hands
- Compression arm sleeves with adjustable bands, with or without hands.

The available compression for the arm sleeves is 23-32 mmHg - Class 2. Available three sizes from 1-S to 3-L in Beige color.




Compression	BASIC ARM	ARM WITH HOOKS	ARM WITH BANDS	Sizes*	Hands	Colors	Composition
CLASS 2 - K2 23-32 mmHg	✓	✓	✓	From 1 to 3	 Without hands  With hands	Beige	70% PA 30% EA



OGUK

GRADUATED COMPRESSION





## BASIC line

Classical line of graduated compression stockings made using high quality double covered elastic fibers, to ensure comfort and durability. Reinforced and anatomic heel for a perfect fit. Four compression degrees available, from 8-11 mmHg (extra light) up to 22-27 mmHg (strong).



COMPRESSION	CODE	ARTICLE		SIZES					COLOURS			COMPOSITION
<b>40 DEN 8-11 mmHg</b>	430	AT	Tights	1	2	3	4XL	5XL	02 Black	36 Skin	21 Soleil	88% PA - 12% EA
<b>70 DEN 12-17 mmHg</b>	750	AD	Knee High Socks	1	2	3	4	5	02 Black	36 Skin	21 Soleil	83% PA - 17% EA
	700	AD	Knee High Socks without Preformed Heel		2	3	4	5	02 Black	36 Skin	21 Soleil	83% PA - 17% EA
	770	AGH	Stay-Up Stockings	1	2	3	4	5	02 Black	36 Skin	21 Soleil	83% PA - 17% EA
	780	AT	Tights	1	2	3	4	5	02 Black	36 Skin	21 Soleil	83% PA - 17% EA
	780X	AT	Tights		2XL	3XL	4XL	5XL	02 Black	36 Skin	21 Soleil	83% PA - 17% EA
	790	ATM	Maternity Tights with Adjustable elastic waist	1	2	3	4	5	02 Black	36 Skin	21 Soleil	83% PA - 17% EA
<b>140 DEN 18-22 mmHg</b>	850	AD	Knee High Socks	1	2	3	4	5	02 Black	36 Skin	21 Soleil	75% PA - 25% EA
	850A	AD	Knee High Socks <b>Open toe</b> <b>NEW</b>	1	2	3	4	5	02 Black	36 Skin		75% PA - 25% EA
	870	AGH	Stay-Up Stockings	1	2	3	4	5	02 Black	36 Skin	21 Soleil	75% PA - 25% EA
	870A	AGH	Stay-Up Stockings <b>Open toe</b> <b>NEW</b>	1	2	3	4	5	02 Black	36 Skin		75% PA - 25% EA
	880	AT	Tights	1	2	3	4	5	02 Black	36 Skin	21 Soleil	75% PA - 25% EA
	880X	AT	Tights		2XL	3XL	4XL	5XL	02 Black	36 Skin	21 Soleil	75% PA - 25% EA
	880A	AT	Tights <b>Open toe</b> <b>NEW</b>	1	2	3	4	5	02 Black	36 Skin		75% PA - 25% EA
	890	ATM	Maternity Tights with Adjustable elastic waist	1	2	3	4	5	02 Black	36 Skin	21 Soleil	75% PA - 25% EA
<b>140 DEN PLUS 22-27 mmHg</b>	850P	AD	Knee High Socks	1	2	3	4	5	02 Black	36 Skin	21 Soleil	75% PA - 25% EA
	870P	AGH	Stay-Up Stockings	1	2	3	4	5	02 Black	36 Skin	21 Soleil	75% PA - 25% EA
	880P	AT	Tights	1	2	3	4	5	02 Black	36 Skin	21 Soleil	75% PA - 25% EA
<b>280 DEN 22-27 mmHg</b>	950	AD	Knee High Socks	1	2	3	4	5	02 Black	36 Skin	21 Soleil	76% PA - 24% EA
	950A	AD	Knee High Socks <b>Open toe</b>	1	2	3	4	5	02 Black	36 Skin		76% PA - 24% EA
	960	AG	Stockings	1	2	3	4	5	02 Black	36 Skin	21 Soleil	76% PA - 24% EA
	960A	AG	Stockings <b>Open toe</b>	1	2	3	4	5	02 Black	36 Skin		76% PA - 24% EA
	970	AGH	DOT Stay-Up Stockings	1	2	3	4	5	02 Black	36 Skin	21 Soleil	76% PA - 24% EA
	970A	AGH	DOT Stay-Up Stockings <b>Open toe</b>	1	2	3	4	5	02 Black	36 Skin		76% PA - 24% EA
	980	AT	Tights	1	2	3	4	5	02 Black	36 Skin	21 Soleil	76% PA - 24% EA
	980X	AT	Tights		2XL	3XL	4XL	5XL	02 Black	36 Skin	21 Soleil	76% PA - 24% EA
	980A	AT	Tights <b>Open toe</b>	1	2	3	4	5	02 Black	36 Skin		76% PA - 24% EA

**Available with the following compression ratings at the ankle**

- 40 DEN 8-11 mmHg** Extra Light, designed for individuals suffering from leg fatigue.
- 70 DEN 12-17 mmHg** Light, designed for individuals with mild circulatory problems, varicose veins and for pregnant women.
- 140 DEN 18-22 mmHg** Medium, designed for individuals with chronic light circulatory problems, varicose veins, swelling and for pregnant women.
- 140 DEN PLUS/280 DEN 22-27 mmHg** Strong, designed for serious circulatory problems, varicose veins, post surgery, oedema, phlebitis and for the prevention of thrombosis.


**Recommendations**

It is recommended to put on the garments in the morning, before getting out of bed. For proper use, please refer to the diagram on the bottom of the box. For proper effectiveness, use the size chart to select the correct size.

**Contraindications**

Individuals with severe arterial circulation disorders, recent and deep vein thrombosis, peripheral artery occlusion and dermatosis with secretions are advised to remove the garment before going to sleep.

SIZE CHART		1	2	3	4	5
Height	cm	154 - 162	158 - 166	162 - 174	166 - 182	174 - 190
Weight	Kg	43 - 53	48 - 63	58 - 73	63 - 83	68 - 93
Ankle girth	cm	18 - 20.5	20.5 - 23	23 - 25.5	25.5 - 28	28 - 30.5
			<b>2XL</b>	<b>3XL</b>	<b>4XL</b>	<b>5XL</b>
Weight	Kg		58 - 68	63 - 78	68 - 88	73 - 98



NEW

# Essential+

*In everything You Do*

## Reinvent the Compression Socks

Essential+ is a revolutionary, high technological content stocking. A **double fabric layer** provides adds elasticity for **improved fit** and enhanced comfort for everyday wear.

The inner **cotton** and **cashmere** lining guarantees a remarkably soft sensation on skin, whereas the outer microfibre layer provides a pleasant sensation of elegance and resistance to the touch.

Product comfort is immediately perceptible during wear: thanks to the **absence of the stretch trim edge** around the knee, **no marks are left on skin** and the fit encourages better circulation.

## OGUK compression stockings



MEDICAL  
DEVICE



DOUBLE  
LAYER



NO LACE  
EFFECT



SOFT  
TOUCH

COMPRESSION	CODE	ARTICLE	SIZES					COLOURS	COMPOSITION
23-32 mmHg - K2	922	AD Knee High Socks	2S	3L	4XL	5XL	6XXL	Black	43% PA - 36% CO - 13% EA - 4% WS - 4% ME

SIZE CHART	2 - S	3 - M	4 - L	5 - XL	6 - XXL
SHOE SIZE	36 - 39	39 - 41	41 - 43	43 - 45	45 - 48
FOOT LENGTH	cm 23 - 25	25 - 26,4	26,4 - 27,7	27,7 - 29	29 - 31
ANKLE GIRTH	cm 18 - 22	20 - 24	22 - 26	24 - 28	26 - 30



CE Dispositivo Medico  
Medical Device  
MADE IN ITALY

OEKO-TEX®  
CONFIDENCE IN TEXTILES  
STANDARD 100  
14CX00025 CENTROCOT  
Tested for harmful substances.  
www.oeko-tex.com/standard100


# COTTON SOCKS

Unisex Cotton knee socks with graduated compression. Fine Cotton at skin contact, soft and comfortable elastic knee hem. Anatomical and reinforced heel and toe for maximum comfort, proper fit and long durability. Two compression degrees available, 18-22 mmHg (medium) and 22-27 mmHg (strong).



COMPRESSION	18-22 mmHg						22-27 mmHg					
CODE	820						920					
ARTICLE	Knee High Socks						Knee High Socks					
COMPOSITION	57% CO - 32% PA - 11% EA						53% CO - 31% PA - 16% EA					
SIZES	1	2	3	4	5	6	1	2	3	4	5	6
COLOURS	01 White		02 Black		01 White		02 Black					
	03 Blue		05 Smoke		03 Blue		05 Smoke					
	15 Grey		36 Skin				15 Grey					

SIZE CHART	1 - XS	2 - S	3 - M	4 - L	5 - XL	6 - XXL
Shoe size	35 - 37	37 - 39	39 - 41	41 - 43	43 - 45	45 - 47
Foot size cm	21.7-23.7	23 - 25	24.4-26.4	25.7-27.7	27 - 29	28.4-30.4
Ankle girth cm	18 - 20	20 - 22	22 - 24	24 - 26	26 - 28	28 - 30



NEW

**NEW 815** Fresh cotton suitable in summer and in the warmer months. Fine Cotton at skin contact, soft and comfortable elastic knee hem.




COMPRESSION	18-22 mmHg					
CODE	815 - FRESH SOCKS					
ARTICLE	Knee High Socks					
COMPOSITION	60%CO-28%PA-12%EA					
SIZES	1	2	3	4	5	6
COLOURS	02 Black		03 Blue			
	05 Smoke					



# COTTON SOCKS EXTRASOFT

Unisex Milk fiber sock with graduated compression. Viscose and Milk fiber make the skin soft and silky and maintains the foot stays fresh and dry. Reinforced and anatomical heel and toe. Product is extremely soft and comfortable. Available in compression 18-22 mmHg (medium).



COMPRESSION	<b>18-22 mmHg</b>					
CODE	820M					
ARTICLE	Knee High Socks					
COMPOSITION	53%VI - 23% PR - 14% PA - 10% EA					
SIZES	1	2	3	4	5	6
COLOURS	01 White			02 Black		
	03 Blue			05 Smoke		


SIZE CHART	1 - XS	2 - S	3 - M	4 - L	5 - XL	6 - XXL
Shoe size	35 - 37	37 - 39	39 - 41	41 - 43	43 - 45	45 - 47
Foot size	cm 21.7-23.7	23 - 25	24.4-26.4	25.7-27.7	27 - 29	28.4-30.4
Ankle girth	cm 18 - 20	20 - 22	22 - 24	24 - 26	26 - 28	28 - 30






# COTTON SOCKS

Cotton and Silver fiber Mens knee sock with graduated compression. Silver fiber at skin contact provides an antibacterial function, and soft and comfortable knee elastic hem. Anatomical and reinforced heel and toe for maximum comfort, proper fit and long durability. Two compression degrees available, 18-22 mmHg (medium) and 22-27 mmHg (strong).



COMPRESSION	18-22 mmHg	22-27 mmHg
CODE	830	930
ARTICLE	Knee High Socks	Knee High Socks
COMPOSITION	60% CO - 22% PA - 10% EA - 8% ME	60% CO - 20% PA - 15% EA - 5% ME
SIZES	2 3 4 5 6	2 3 4 5 6
COLOURS	02 Black	02 Black

SIZE CHART	2 - S	3 - M	4 - L	5 - XL	6 - XXL
Shoe size	37 - 39	39 - 41	41 - 43	43 - 45	45 - 47
Foot size cm	23 - 25	24.4 - 26.4	25.7 - 27.7	27 - 29	28.4 - 30.4
Ankle girth cm	20 - 22	22 - 24	24 - 26	26 - 28	28 - 30



OGUK

DIABETIC

OGUK

## DIABETIC AND SENSITIVE FEET SOCKS

Unisex Cotton socks with Silver fiber for diabetic and sensitive feet, without compression and seamless. Silver, thanks to its beneficial natural properties, provides antibacterial and thermoregulating actions, keeping the feet in a protected and healthy environment. Reinforced and anatomical toe and heel design for perfect fit and long durability. Available in three different lengths and also available with soft sole.



x-static®



# DIABETIC AND SENSITIVE FEET SOCKS

Unisex Cotton socks with Crabyon fiber for diabetic and sensitive feet. With no compression and seamless design. Crabyon fiber, along with its natural properties of Chitin, provides protection and a regenerating action which maintains the foot stays in a protected and healthy environment. Reinforced and anatomical toe and heel design, offers a perfect fit and long durability. Available in three different lengths.



OGUK

## DIABETIC AND SENSITIVE FEET TOE SOCKS

Unisex Cotton toe socks with Silver fiber for diabetic and sensitive feet. No compression and seamless construction. Silver, along with its natural properties, provides antibacterial and thermoregulating actions, keeping the feet protected and in a healthy environment. The construction with individual toe sleeves, keeps the area clean and healthy. Available in 2 different lengths.



x-static®



OGUK

ORTHOPAEDIC RANGE

# OGUK

## ORTHOPAEDIC WRISTBAND AND MITTEN

These articles are normally recommended for all cases where the immobilization of the hand articulation is required, after operation, in case of arthritis, after plaster bandages, during the rehabilitation period and for sport inflammation.



Art. 40100



Art. 40150

### Artt. 40100 • 40150

Size	
1^	12/ 17
2^	18/ 23
3^	24/ 30



### Art. 40100

Orthopaedic wristband with Velcro regulation in col. Beige

**Composition:**  
54% Polyester (PL), 20% Nylon (PA), 20% Cotton (CO), 6% Elastane

### Art. 40150

Orthopaedic wristband with Velcro regulation in col. White

**Composition:**  
54% Polyester (PL), 20% Nylon (PA), 20% Cotton (CO), 6% Elastane

Art. 40200



### Art. 40200

Size	
1^	11/ 15
2^	16/ 22
3^	23/ 28
4^	29/ 34



### Art. 40200

Elastic orthopaedic mitten.

**Composition:**  
60% Cotton (CO), 32% Nylon (PA), 8% Elastane (EA)



Art. M1800



Art. M1900

### Artt. M1800 • M1900

Size	
S	15-17
M	17-19
L	19-21
XL	21-23
XXL	23-25



### Art. M1800

Elastic orthopaedic mitten with support aluminium insert.

**Composition:**  
41% Polyester (PL), 32% Elastane (EA), 27% Metal (ME)

### Art. M1900

Elastic orthopaedic mitten with support aluminium insert and closing with Velcro.

**Composition:**  
41% Polyester (PL), 32% Elastane (EA), 27% Metal (ME)

# OGUK

**LINE OF ELASTIC SUPPORT TUBULAR  
BANDAGES FOR KNEE, ANKLE, CALF  
AND THIGH**

Articles of the strong (class 1) and extra strong (class 2) classes of compression. Articles are indicated in all cases where fixation and the compression are required: after breaks, tendon and muscles tension, during rehabilitation period after plaster bandages with the immobilization of articulation, after operations...

The articles have in the inside surface the cotton, which makes them hypo-allergenic and ensures soft contact with the skin. Because of the special type of production, the articles are very thin, but also durable.





Art. 50100



Art. 50200



Art. 50300

**Art. 50100**  
Elastic ankle support band with compression class II.  
**Composition:**  
57% Nylon (PA), 23% Elastane (EA), 20% Cotton (CO)

**Art. 50200**  
Elastic knee support band with compression class II.  
**Composition:**  
57% Nylon (PA), 23% Elastane (EA), 20% Cotton (CO)


**Art. 50300**  
Elastic calf support band with compression class II.  
**Composition:**  
57% Nylon (PA), 23% Elastane (EA), 20% Cotton (CO)

**Art. 50100**




Size	cm
1^	19/ 21
2^	21/ 23
3^	23/ 26
4^	24/ 29
5^	29/ 33

**Art. 50200**



Size	cm
1^	29/ 33
2^	33/ 37
3^	37/ 41
4^	41/ 44
5^	44/ 49

**Art. 50300**



Size	cm
1^	24/ 28
2^	28/ 32
3^	32/ 36
4^	36/ 40
5^	40/ 44

**LINE OF ELASTIC SUPPORT TUBULAR BANDAGES FOR KNEE, ANKLE, CALF AND THIGH**




Art. 70100



Art. 70300

**Artt. 70100 70300**



Size	cm
1^	19/ 21
2^	21/ 23
3^	23/ 26
4^	24/ 29
5^	29/ 33

**Art. 70100**  
Elastic knee socks support band with compression class I.  
**Composition:**  
60% Cotton (CO), 23% Nylon (PA), 17% Elastane (EA)

**Art. 70300**  
Elastic ankle support band with compression class I.  
**Composition:**  
60% Cotton (CO), 23% Nylon (PA), 17% Elastane (EA)



Art. 70200

Art. 70200	
Size	cm
1^	29/ 33
2^	33/ 37
3^	37/ 41
4^	41/ 44
5^	44/ 49



Art. 70400

Art. 70400	
Size	cm
1^	24/ 28
2^	28/ 32
3^	32/ 36
4^	36/ 40
5^	40/ 44



Art. 70500

Art. 70500		
Size	cm	cm
1^	36-37,5	44/ 45,5
2^	39-40	46,5/ 48
3^	41,5-43	50/ 51,5
4^	44,5-45,5	53/ 55
5^	47-48,5	57/ 59

**Art. 70200**  
Elastic knee support band with compression class I.  
**Composition:**  
60% Cotton (CO), 23% Nylon (PA), 17% Elastane (EA)

**Art. 70400**  
Elastic calf support band with compression class I.  
**Composition:**  
60% Cotton (CO), 23% Nylon (PA), 17% Elastane (EA)

**Art. 70500**  
Elastic thigh support band with compression class I.  
**Composition:**  
60% Cotton (CO), 23% Nylon (PA), 17% Elastane (EA)



Art. C4400

Art. C4400	
Size	cm
S	22-25
M	25-27
L	27-31
XL	31-34
XXL	34-37

**Art. C4400**  
Elastic ankle support band with lateral silicone insert.  
**Composition:**  
55% Nylon (PA), 45% Elastane (EA)

**LINE OF ELASTIC SUPPORT TUBULAR BANDAGES FOR KNEE, ANKLE, CALF AND THIGH**

# OGUK

## KNEE AND ARMS SUPPORT

They are used in case of pains in the knee joint, instability, arthritis, after operations at the knee articulation and in all cases where is necessary the stabilization of the knee articulation.  
Line of orthopaedic knee-bands and arm supports.



Art.  
G2700

Art.  
G2800

Artt. G2700 G2800	
Size	cm
S	34-36
M	36-38
L	38-41
XL	41-44
XXL	44-47

### Art. G2700

Knee guard with two lateral flexible strips for better stability. Neoprene 5 mm.  
**Composition:**  
80% Neoprene, 15% Cotton (CO), 5% Nylon (PA)

### Art. G2800

Knee guard with Velcro fastener. Neoprene 5 mm.  
**Composition:**  
80% Neoprene, 15% Cotton (CO), 5% Nylon (PA)

Art.  
B3000



Art. B3000	
Size	cm
S	30-35
M	35-40
L	40-45
XL	45-50
XXL	>50

### Art. B3000

Arm support in order to immobilize the shoulder and to avoid abduction and rotation.  
**Composition:**  
Stratified, cotton with foam rubber padding.



# OGUK

## Contact Us

---

T : 0116 260 9117

E : [info@oguk.co.uk](mailto:info@oguk.co.uk)

W : [www.oguk.co.uk](http://www.oguk.co.uk)

OGUK Ltd, Unit 1, Wanlip Road, Syston, Leicester LE7 1PD